



Flash Receiver BE1442/1443

Flash Receiver BE1441

Settings

No adjustments are required for normal use. The relevant descriptions are provided below, if you wish to change a setting for some reason.

Radio key

In order to use several Bellman Visit Systems close to one another without interference, different Radio Keys can be set on the different systems. All Bellman Visit System units are supplied from the factory tuned to the same Radio Key, channel 0. This means that all Radio Key Switches on the transmitters are set to the OFF position.

A Bellman Visit transmitter is also required to alter the radio key on the Bellman Visit Flash Receiver. Proceed as follows to change the Radio Key:

- Set a Bellman Visit Transmitter to the desired Radio Key by altering its Radio Key Switch to the desired Radio Key. See the transmitter user manual for further information.
- Hold down the Test Button (9) on the Bellman Visit Flash Receiver until the green  and yellow  LEDs blink alternately. The Bellman Visit Flash Receiver will now be in programming mode for about 30 seconds.
- Press the Bellman Visit transmitter's test button.
- The LEDs     will now blink five times in quick succession to indicate a successful change of radio key.
- After changing the radio key, the Bellman Visit Flash Receiver will automatically return to normal mode.

Please note: all Bellman Visit products within a system must be tuned to the same Radio Key in order to operate as a group.

Signal pattern

A Signal Pattern is the name for the way in which a receiver in the Bellman Visit System indicates activation. Changing the transmitters' Signal Switch changes the Signal Pattern which the receivers display when the transmitter is activated.

The following signal patterns are available for the Bellman Visit System: (Please note that the Flash Receiver does not emit any sound):

868/315MHz Europe, Asia


Type	LED pattern	Sound	Vibration	Flash
Green 1	Green is constantly lit	1 x ding dong, low-frequency tone	Separate	Yes
Green 2	Green blinks in sequences of two	2 x ding dong, low-frequency tone	Separate	Yes
Green 3	Green blinks in sequences of three	1 x ding dong, high-frequency tone	Separate	Yes
Green 4	Green blinks constantly	2 x ding dong, high-frequency tone	Separate	Yes
Yellow 1	Yellow is constantly lit	1 x ring, low-frequency tone	Short	Yes
Yellow 2	Yellow blinks in sequences of two	2 x ring ring, low-frequency tone	Short	Yes
Yellow 3	Yellow blinks in sequences of three	1 x ring, low-frequency tone	Short	Yes
Yellow 4	Yellow blinks constantly	2 x ring ring, high-frequency tone	Short	Yes
Orange 1	Orange is constantly lit	Baby	Rapid	Yes
Orange 2	Orange blinks in sequences of two	Baby	Rapid	Yes
Orange 3	Orange blinks in sequences of three	Baby	Rapid	Yes
Orange 4	Orange blinks constantly	Baby	Rapid	Yes
VMA	Red and Orange constantly blink alternately	VMA constant	Long	Yes
Fire alarm	Red blinks constantly	Fire alarm constant	Long	Yes

433MHz North America, Australia

Type	LED pattern	Sound	Vibration	Flash
Green 1	Green is constantly lit	1 x Ding dong, Bass tone	Long	Yes
Green 2	Green is constantly lit	1 x Ring, Treble tone	Short	Yes
Green 3	Green is constantly lit	1 x Ding dong, Bass tone	Long	Yes
Yellow 1	Yellow is constantly lit	1 x Ding dong, Treble tone	Long	Yes
Yellow 2	Yellow is constantly lit	1 x Ring, Bass tone	Short	Yes
Yellow 3	Yellow is constantly lit	1 x Ding dong, Treble tone	Long	Yes
Orange 1	Orange is constantly lit	2 x Ding dong, Bass tone	Long	Yes
Orange 2	Orange is constantly lit	2 x Ring, Bass tone	Short	Yes
Red 1	Red is constantly lit	2 x Ding dong, Treble tone	Long	Yes
Red 2	Red is constantly lit	2 x Ring, Treble tone	Short	Yes
Red 3	Red is constantly lit	2 x Ding dong, Treble tone	Long	Yes
VMA	Red and Orange constantly blink alternately	VMA	Constant	Yes
Fire alarm	Red blinks constantly	Fire	Constant	Yes


Battery backup

BE1442

Getting started

- Remove the battery slip on the bottom of the unit before use.
- Please charge the unit for 24 hours to ensure full backup time.


Replacing batteries (Done by the operator)

- Remove the screw on the bottom of the unit to exchange the batteries.

CAUTION

RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.

DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.

- Use only 4x1.2VNiMH 600AAA/900AAA Lexel rechargeable batteries to ensure full backup time.

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Advanced programming

Advanced programming provides additional options for those who wish to make special modifications to the Bellman Visit Flash Receiver.

The idea is that it should be possible to select a completely unique signal pattern which is linked to activation from a specific input on a special Bellman Visit transmitter. The function works regardless of the radio key settings on the units that are programmed. Please note that, for safety reasons, the function does not work with the BE1480 Bellman Visit Smoke Alarm. By using advanced programming of the Bellman Visit Flash Receiver, it can be adjusted so that its signal pattern corresponds exactly to what is required. In other words an entirely individual signal pattern can be programmed, such as displaying an orange permanently lit LED and a constant vibration.

In order to adjust the setting, the Bellman Visit Transmitter to which the Bellman Visit Flash Receiver should be adapted must be available. The transmitter must also be connected so that it can be activated in the way in which it is intended to be used.

Proceed as follows:

- 1 Hold down the Test Button (9) on the Bellman Visit Flash Receiver until the green  and yellow  LEDs blink alternately. The Bellman Visit Flash Receiver will now be in programming mode for about 30 seconds.
- 2 Hold down the Test Button (9) at the same time as the relevant Bellman Visit transmitter is activated in precisely the way in which it is intended to be used. Note that all inputs are individual. It is therefore not possible to use the Test Button on a Bellman Visit Telephone Transmitter whose telephone input will indicate the relevant pattern.
- 3 Scroll through the different LED options by a short press on the Test Button (9). Select the relevant indication by holding down the Test Button (9) until the LED  goes out and starts to shine with a constant green light again.
- 4 Scroll through the different vibration options by a short press on the Test Button (9). Select the relevant indication by holding down the Test Button (9) until the LED  goes out and starts to shine with a constant green light again.
- 5 The Bellman Visit Flash Receiver will now show the indication method programmed. End the display with a short press of the Test Button (9).
- 6 After a short while, the Bellman Visit Flash Receiver will automatically return to normal mode. This function is essential where a transmitter has to work in a special way with regard to a specific receiver.

Resetting advanced programming

It is quite easy to reset the Bellman Visit Flash Receiver if it needs to be reset after it has been programmed using advanced programming.

- 1 Hold down the Test Button (9) on the Bellman Visit Flash Receiver until the green  and yellow  LEDs blink alternately. The Bellman Visit Flash Receiver will now be in programming mode for about 30 seconds.
- 2 Press the Test Button (9) three times in quick succession.
- 3 All LEDs     remain constantly on for a few seconds.
- 4 All the advanced programming has now been deleted and the Bellman Visit Flash Receiver will automatically return to normal mode.

Testing

It is easy to test the BE1441/1442/1443 Bellman Visit Flash Receiver. If the Flash Receiver does not work as described below, you can check further below under Appendix/Further information/Troubleshooting/Troubleshooting guide.

How to test

A transmitter in the Bellman Visit System which is set to the same Radio Key as the Bellman Flash Receiver is required to test the flashing light and vibration, if the BE1270 Bellman Bed-shaker (accessory) is connected, and the radio reception on the BE1441/1442/1443 Bellman Visit Flash Receiver.

- Press the transmitter test button.
- The Bellman Visit Flash Receiver will give the following indications:
 - o Flashing with the flash light (2).
 - o The LEDs     which the transmitter has been set to indicate with or which have been programmed into the Flash Receiver with Advanced Programming will be switched on.
 - o If the BE1270 Bellman Bed-shaker (accessory) is connected, it will vibrate as the transmitter has been set to indicate or in the way the Flash Receiver has been programmed with Advanced Programming.

To test the analogue telephone input:

- Connect the Bellman Visit Flash Receiver to an analogue telephone socket via the Telephone Input . Use the BE9105 Telephone Flex (accessory) and an adapter plug (accessory).
- Ring the telephone number. The Bellman Visit Flash Receiver will then start flashing, the yellow LED  will light up and, if a BE1270 Bellman Bed-shaker (accessory) is connected, it will vibrate.

A short press on the Test Button (9) activates BE1441/1442/1443 so that it repeats its last indication.

Troubleshooting

You can carry out a number of checks yourself before sending a product for repair.

Problem	Solution
Nothing happens.	<ul style="list-style-type: none"> Check that the power supply unit is connected correctly. The LED ⚡ should light up green. Check that there is current in the wall socket.
The receiver is not activated.	<ul style="list-style-type: none"> Check the battery in the transmitters. Check that the receiver is not placed too far away from the transmitters by moving it closer to the transmitters. Check that the receiver is set to the correct radio key. For further information see Function/Radio key.
The receiver signals when no transmitter is activated.	<ul style="list-style-type: none"> Change the Radio Key on all units in the system. For further information see Function/Radio key.

Technical information

Power supply

Mains power:

7.5 V DC / 1500 mA.

This product is also designed for IT power distribution system with phase to phase voltage of up to 230V. The socket-outlet shall be installed near the equipment and shall be easily accessible.

Europe: BE9092 (Something High Electric P12-075150EU)

Great Britain: BE9217 (Something High Electric P12-075150UK)

US/CAN: Something High Electric P12-075150US

Australia/NZ: Something High Electric P12-075150AS

Battery back-up

for BE1441: not applicable

for BE1442: 4x1.2V Lexel 600/900mAh

for BE1443: 4x1.2V Lexel 900mAh

- Please dispose used batteries at municipal collecting points or dispense them to local stores free of charge.

Radio function

Power consumption:	Active: 1500 mA
Radio frequency:	315/433.92/868.3 MHz depending on region
Number of Radio Keys:	64 Radio Keys as standard. Special software can be used to increase these to 256 Radio Keys in increments of 64 per software purchase. Contact the nearest supplier for further information.
Coverage:	The normal coverage between a transmitter and receiver in the Bellman Visit System is up to 200 metres (868MHz)/ 80m (433MHz)/ 50m (315MHz) with clear line of sight. Coverage is reduced if walls and large objects screen off the signals. Any thick walls constructed of reinforced concrete will greatly affect coverage. The system may also be affected by radio transmitters such as TV transmitters, computers, mobile phones, etc. This means that a unit may work perfectly in one part of the room but not at all in another.

Additional information

Programming pin (16):	For service personal only (Output: 3.3V/max. 150mA) CAUTION: Do not use, unit will be damaged when used.
-----------------------	--

Activation

Radio:	Bellman Visit system
Via analogue telephone network:	26 - 120 V RMS, 15 - 100 Hz.
BE1441/BE1442:	Via mobile phone detector  : see accessories (Output: 3.3V/max 150mA)

WARNING: The telephone port of the Flasher Receiver must be connected using 26 AWG or greater telecommunications wire.

Output signals

Output signals	Built-in flash light signal >10 Candela Warning! Flashes can cause epileptic attacks.
Vibrator power	2.0 – 4.0 VDC (max. 500mA)
For indoor use only	Temperature range 0°C – 40°C
Dimensions Ø x H:	70 x 140 mm
Weight:	320 g
Colour:	White

Optional accessories (sold seperately)

Wall Mount Bracket BE9075

Bellman Bed-shaker BE1270

Telephone Flex BE9105

Mobile Detector BE9250

Adapter plug for the appropriate country

Hereby, Bellman & Symfon, declares that this wrist receiver and charger is in compliance with the essential requirements and other relevant provisions of directives: R&TTE 1999/5/EC, LVD 1973/23/EC, EMC 1989/336/EC and MDD 1993/42/EEC.

For use in all EU countries

Correct use and disposal of batteries. Replace only with the same or equivalent type recommended by the manufacturer. Please dispose of old, defective batteries in an environmentally friendly manner in accordance with the relevant legislation. Dispose used batteries at municipal collecting points or dispense them to local stores free of charge.

Contact us

For a complete Declaration of Conformity please contact info@bellman.com or call +46 31 68 28 20

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Bellman & Symfon Europe AB
Bellman & Symfon Asia Ltd.