Wisbi®

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
SAFETY INFORMATION

Please carefully read through these instructions before using the device! TEXIBLE GmbH assumes no liability for material damages and personal injuries that arise due to non-observance of the operating instructions and especially the safety information!

TEXIBLE Wisbi® is not a medical device and must not be used to diagnose or treat illnesses. Ask your doctor for prevention, diagnosis and treatment of illnesses.

The TEXIBLE Wisbi transmitter must be placed under or next to the bed in a way that there is no danger for the patient through the cable.

Note for care workers

TEXIBLE Wisbi is a support aid for care workers and does not serve as a replacement. TEXIBLE Wisbi must be checked for functionality before each use.

Maintenance and troubleshooting

TEXIBLE Wisbi is free of maintenance except when the battery needs replacement. Only tighten the screws on the battery compartment slightly.

The following steps should be considered if the function is restricted:

- Is the bed pad correctly connected and the transmitter switched on?
- Is the TEXIBLE Wisbi transmitter connected with the receiver?
- Is the transmitter battery functioning?
- Is the bed pad free of defects?

If the mentioned points are fulfilled, but the function is still limited, please contact the technical staff.

Note: The TEXIBLE Wisbi transmitter and the receiver form an encrypted paging system. Both should be clearly marked that they don't get mixed up. This ensures that the call is only triggered for the room in question.

Make sure that the radio link is unobstructed. Do not mount in distribution boxes, metal enclosures or close to large metal objects.

1. PACKAGE CONTENTS

The package contents:

- 1x TEXIBLE Wisbi® – Bed Pad
- 1x TEXIBLE Wisbi® Plus – Transmitter (compatibel with ELDAT EasyWave® Receiver)

2. SCOPE OF APPLICATION AND FUNCTION

TEXIBLE Wisbi is the innovative nursing care product that stands out for its ability to automatically signal an alarm when moisture, body fluids or a bed exit are detected.

The sensor textile is used in a wide variety of nursing applications, such as senior and nursing homes, assisted living and family care. The necessary basic care, as well as the incontinence care and the supervision of demented and fall-endangered residents represents an increased care effort for all persons involved.

TEXIBLE Wisbi offers the right support here. The system allows an optimal skin climate and minimizes the risk of running away.

3. PUTTING TEXIBLE WISBI INTO OPERATION

The TEXIBLE Wisbi transmitter is equipped with an 868MHz EasyWave® Radio transmitter and forms a call system with an ELDAT receiver. Every 2 seconds a measurement is made. Upon detection of a wet bed pad or a bed exit, an alarm is sent via radio.

- 1. INSERT BATTERIES:** Insert three mignon batteries (AA) into the battery compartment. When the batteries are inserted correctly, the white LED flashes briefly.

Switch the TEXIBLE Wisbi transmitter on through pressing the button  once. The white LED flashes.

- 2. CONNECTION BED PAD:** Place the bed pad as the top layer on the sheet. White cotton side up, blue side down. Then fasten the TEXIBLE Wisbi transmitter using the push buttons with the bed pad. (See image below)

3. CONNECT RECEIVER WITH TRANSMITTER:

Put the Eldat Receiver into operation in accordance with the accompanying operating instructions and activate the „teach-in mode“ of the receiver.

Pair the TEXIBLE Wisbi Plus Transmitter with the ELDAT receiver within 20 seconds by pressing the button  for 3 seconds.

4. ADDRESS CHANGE OF THE TRANSMITTER:


If several EasyWave radio transmitters are in use, it may be the case that the factory-set addressing of the TEXIBLE Wisbi transmitter is identical to other devices. In this case, the radio address of the TEXIBLE Wisbi transmitter must be changed. To do so, open the 4 screws on the rear support cover. The address can then be changed in the viewing window via the DIP switch (Number 1 - 4).

You can now close the support cover once again and couple the TEXIBLE Wisbi transmitter as described in Section 3: „Connect Receiver with Transmitter“.

5. ALARM SETTINGS OF TEXIBLE WISBI TRANSMITTER:

The alarm settings (occupancy detection ON/OFF, delays) of the transmitter can be customized to your Needs:

To do so, open the 4 screws on the back Support coverage. The address can then be changed in the inspection window via the DIP switch (Number 5 and 6).

Alarm WET and OCCUPIED immediately

Alarm WET immediately and OCCUPIED deactivated

Alarm WET immediately and OCCUPIED 5 minutes delayed

Alarm WET immediately and OCCUPIED 15 minutes delayed

Then close the support cover again.

Functional test

TEXIBLE Wisbi must be checked for proper functioning before each use. Switch on the TEXIBLE Wisbi transmitter. Press the button  for 3 seconds until the orange LED flashes. This will lead to a functional test of the electronics and sends an alarm.

Switching TEXIBLE Wisbi off

Once TEXIBLE Wisbi is switched on it is constantly in standby mode and measures every 2 seconds.

If the transmitter is not used, it can be switched off by pressing the button  for 3 seconds to preserve the batteries. By pressing the button  again once, the transmitter is switched on again.

Changing batteries of TEXIBLE Wisbi

Open the rear cover of the battery compartment by loosening the screw. Insert three mignon batteries (AA) into the battery compartment as shown and close the cover again.

4. OPERATING CONCEPT

Function	Button	LED Display
Switching TEXIBLE Wisbi Transmitter ON	 Press 1 x	flashing white
Switching TEXIBLE Wisbi Transmitter OFF	 Press 3 sec	glowing red
Functional test	 Press 3 sec	flashing orange
Alarm		flashing orange
Battery low		red & orange flash alternately
Error		glowing red

5. DISPOSAL INFORMATION

Dispose the TEXIBLE Wisbi transmitter with accompanying batteries at a designated collection point.

The TEXIBLE Wisbi bed pad can be disposed in household waste. Please contact your waste disposal company or the representative for more information about nearby collection points.

Dispose of the packaging material in the collection containers for cardboard, paper and plastics.

6. WARRANTY

We rectify any device defects due to material or manufacturing faults by repair or replacement free of charge within the statutory warranty period. The warranty expires in the case of intervention by third parties.

7. EC - DECLARATION OF CONFORMITY

Manufacturer:
TEXIBLE GmbH, Schwefelbadstrasse 2, 6845 Hohenems, Austria

The following guidelines have been complied with:
EMC 2014/30/EU, RED 2014/53 EU, LVD 2014/35/EU, RoHS 2011/65/EU


The declaration of conformity can be requested from the manufacturer.

8. CUSTOMER SERVICE

Please contact the manufacturer or your dealer if malfunctions still occur despite proper handling or if the device has been damaged.

9. TECHNICAL DATA TEXIBLE WISBI

Size	64 x 100cm
Weight	250g
Water absorption	up to 700ml/m²
Trigger quantity	10ml
Minimal body weight (for occupancy detection)	45kg
Material	67% Cotton 30% Polyester 3% Polyurethane
Wash resistance	    
Power supply	3xAA mignon batteries, LR6 -----
Power	1,4mW
Protection class	IP40
Ambient temperature	0 to +55°C
Battery life	1 year depending on application
Radio	Eldat EasyWave® 868,3 MHz
Radio range	Up to 30m in the building, depending on construction conditions

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TEXIBLE GmbH, Schwefelbadstrasse 2, 6845 Hohenems, Autriche