

Boost your hearing aid
Connect a neck loop, and the Mino further amplifies speech and acts as an additional volume and tone control for your hearing aid.



Discover Mino

The Mino is a new speech amplifier that helps you hear better in most everyday situations, whether you wear a hearing aid or just suffer the occasional hearing difficulty.



Sound with CD quality
The Mino uses similar digital technology as modern CD/DVD players and delivers crystal clear speech, free from background noise.



Pocket size and super light
The Mino is as small as a modern cellphone and weighs just 2½ oz. So the Mino can be kept discreetly in your pocket without weighing you down.



Hearing solutions for everyday life

For the past 20 years, Bellman & Symfon of Sweden has worked hard to improve quality of life for people with hearing impairments. Our hearing products comprises solutions that have been thoroughly researched and developed in close collaboration with audiologists, users and leading experts. This guarantees products with great sound quality and unparalleled ease of use.

Distributor:
Seniorentechnik-Martin.de
Martin Kommunikationstechnik
Bockhorster Landweg 30a
33775 Versmold
Tel. 05423/4734-0
Fax. 05423/473420
eMail: info@seniorentechnik-martin.de
Homepage: www.seniorentechnikmartin.de
Shop: www.Seniorentechnik-Martin.de



Head Office
Bellman & Symfon Europe AB
Södra Långebergsgatan 30
421 32 Västra Frölunda
Sweden
Phone +46 31 68 28 20
Fax +46 31 68 28 90
E-mail info@bellman.se
www.bellman.com

China Office
Bellman & Symfon Asia Ltd.
1115-1116, Main Tower, GD Int'l Hotel
339 Huanshi Dong Rd
Guangzhou, China, 510095
Phone +86 20 8318 2950
Fax +86 20 8333 5528
E-mail info.cn@bellman.com
www.bellman.cn

© 2010 Bellman & Symfon AB. Bellman, Bellman Audio and Bellman & Symfon are registered trademarks owned by Bellman & Symfon AB, all rights reserved. All Product designs are patented.

CL_00185H001_EN

Bellmanaudio
FOR EASY COMMUNICATION

Mino Professional communicator



Having difficulty hearing?



Are these situations familiar?


- At home.** Does your partner often need to speak up to get her point across?
- With friends.** Do you have trouble following the conversation when you are out with friends?
- In theaters.** Do you avoid going to the movies and concerts?
- In stores.** Do you have trouble understanding what the clerk says when out shopping?
- Watching TV.** Do you miss words when you watch TV at normal volume?
- With hearing aid.** Do you feel that your hearing aid doesn't fully help in some situations?

Discover the easy way to hear



Single touch action
The Mino is intuitive to use. All functions are accessed with the simple press of a button and there are no complicated menus to navigate.

Long lasting power
The Mino uses a rechargeable lithium-ion battery, just like a modern cell-phone. So you never need to worry about buying new batteries.




Zoom the sound
Mino has both an omni- and a unidirectional microphone (zoom). With the zoom microphone, you can 'point' to the person talking, reducing any background noise.

Bring Mino to the movies
The Mino is equipped with a telecoil. The telecoil picks up sound in venues fitted with hearing loop systems, such as theaters, churches, and concert halls.



Customize your Mino
If you'd like to tailor your solution, we offer a number of accessories, ranging from practical car chargers and travel cases to advanced hearing accessories. It's your decision.

Discover the easy way to hear better with digital technology




At home



Listen at the volume of your choice without others having to raise their voices. Simply place the Mino near the other person and worry no more.

The Mino effectively isolates speech while removing disturbing background sounds. The result is crystal clear speech with no noise. We offer a number of accessories for the Mino. Non-hearing aid users can, for example, choose to listen using headphones, MicroSet™ earphones or stetoclips.

▼ Headphones ▼ Earphones ▼ Stetoclips



With friends



Trouble hearing people in noisy environments? Change to the zoom microphone to reduce background noise and focus on the conversation.




In theaters



Do you avoid movie theaters and concerts? Many public venues, such as theaters, concert halls, and churches, are fitted with hearing loops.

Sound from the loudspeakers in such venues is converted into an electromagnetic field that the Mino picks up via its built-in telecoil. The Mino's telecoil always delivers perfect sound, whether it's placed in your pocket or on a table. Look for a T symbol sign next time you visit your favorite movie theater.



In stores



When you're shopping down town, it's handy to keep the Mino in a pocket so that both your hands are free.

The MicroSet™ is a unique accessory specially developed for such times. It includes a small lapel microphone that you simply attach to your clothes and a discreet earphone that you pop in your ear. The lapel microphone is unidirectional (zoom) and reduces background noise in the store, so you can talk with the clerk uninhibited.


▼ MicroSet™



▼ USB charger*



* Included in the package



Watching TV




When you have difficulty hearing, it's easy to turn the TV volume up to the max, which other family members don't always appreciate.

With the Mino, there's no reason to turn up the TV volume to hear better. Use a cable kit to connect the Mino to the TV and you can enjoy great TV sound in stereo while sparing your family's ears. If you want to comment on the show, simply press the Push-to-Listen button on the Mino to lower the volume and join the conversation.

▼ Cable kit, 16 feet



▼ Car charger



With hearing aid



If your hearing aid isn't quite enough in certain situations, you should try connecting it to your Mino using our discreet neck loop.

The Mino further amplifies the sound and also acts as a remote for easy adjustment of the volume and treble, as well as microphone selection. You can also connect the Mino to the TV to amplify the TV sound in your hearing aid. Do you travel a lot? Then you can store your Mino and accessories in this protective travel case for the Mino.

▼ Neckloop



▼ Travel case

